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High performance liquid chromatography (HPLC)

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Automated High Pressure Gradient System
 Suitable for R&D work, providing high accuracy and high sensitivity for a variety of samples.

٢٢١


SIL-10Axl
 Auto Injector
 provides high injection volume repeatability and low cross contamination
 (*RSD value 0.3%)

CBM-10A
 Communications Bus Module
 *with reliable optical fiber interface (for CLASS-LC workstation. The SCL-10A System Controller is available for other workstation or integrators)

SPD-10A/10AV
 UV-VIS Detector
 provides an exceptionally low noise level and simultaneous dual wavelength detection
 (*noise level: $\pm 0.35 \times 10^{-4}$ AU)

Reservoir Tray
 *accommodates 7 one-liter bottles

CTO-10A/10AC
 Column Oven
 provides precise temperature control of columns, mixers and switching valves.
 (temperature control precision: $\pm 0.1^\circ\text{C}$)

Option Box L
 *accommodates small components and peripherals

Sample Cooler
 *keeps sample vials at low temperature

DGU-10B
 Helium Degassing Unit
 *improves reliability of high sensitivity analysis

LC-10AD
 Solvent Delivery Unit
 *provides low pulsation solvent delivery

CLASS-LC10
 Chromatography Workstation
 *supports GLP/GMP compliance

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
BONDED PHASES FOR HPLC AND THEIR ABBREVIATIONS

Phase	Description	Phase	Description
S1	Silica $-\text{Si}-\text{OH}$ Classic normal phase material. Suitable for separating polar non-ionic organic compounds.	CN	CPS, PCN, Cyano, Cyanopropyl, Nitrile $-\text{Si}-\text{CH}_2\text{CH}_2\text{CH}_2\text{CN}$ Can be employed as either a Reversed Phase or normal phase material. Slightly polar, unique selectivity for polar compounds in both Reversed Phase and normal phase modes. Equilibrates very rapidly, suitable for gradient separations. Useful for many pharmaceutical applications (e.g. tricyclic antidepressants).
C1	TMS, SAS, Trimethyl $-\text{Si}-\text{CH}_3$ Reversed Phase material. Unique selectivity for polar and multifunctional compounds. Least retentive of all alkyl group bonded phases for non-polar solutes.	NH ₂	APS, Amino, Amino Propyl Silyl $-\text{Si}-\text{CH}_2\text{CH}_2\text{CH}_2\text{NH}_2$ Can be employed as Reversed Phase, normal phase or weak anion exchange material. Reversed Phase useful for separating carbohydrates. Normal phase alternative selectivity to silica, not deactivated by small amounts of water. Ion exchange: weak anion exchanger when used with buffers. Separates anionic and organic acids.
C2	RP-2, Dimethyl $-\text{Si}-\text{C}_2\text{H}_5$ Reversed Phase material. Less retentive than C4, C8, or C18. More retentive than C1.	NO ₂	Nitro $-\text{Si}-\text{NO}_2$ Normal phase material. Separates aromatic compounds and compounds with double bonds.
C3	Propyl $-\text{Si}-\text{C}_3\text{H}_7$ Reversed Phase material. Used in Hydrophobic Interaction Chromatography (HIC) of proteins and peptides.	OH	Diol, Glycerol $-\text{Si}-\text{CH}_2\text{CH}_2\text{CH}_2\text{OCH}_2\text{CH}(\text{OH})\text{CH}_2\text{OH}$ Can be employed as either a Reversed Phase or normal phase material. Reversed Phase: used for Gel Filtration Chromatography (GFC) of proteins and peptides. Normal phase: similar selectivity to silica not deactivated by small amounts of water.
C4	Butyl $-\text{Si}-\text{C}_4\text{H}_9$ Reversed Phase material. Useful for ion-pairing chromatography. Offers less retention than C8 and C18 phases for non-polar solutes. When bonded to 300Å silica, it is an ideal phase for analyzing large proteins and hydrophobic peptides.	SAX	SB, Quarternary amine, Strong Base $-\text{Si}-\text{CH}_2\text{CH}_2\text{CH}_2\text{N}^+(\text{CH}_3)_3$ Ion exchange material. Strong anion exchangers (basic) are useful for separating nucleotides, nucleosides, and organic acids.
C6	Hexyl $-\text{Si}-\text{C}_6\text{H}_{13}$ Reversed Phase material. Useful for ion-pairing chromatography. Less retentive than C8 and C18 phases.	SCX	SA, Sulfonic Acid, Strong Acid $-\text{Si}-\text{CH}_2\text{CH}_2\text{CH}_2\text{SO}_3\text{H}$ Ion exchange material. Strong cation exchangers (acidic) are useful for separating organic bases.
C8	MOS, RP-8, LC8, Octyl $-\text{Si}-\text{C}_8\text{H}_{17}$ Reversed Phase material. Similar selectivity to C18 but less retentive. Wide applicability (e.g. pharmaceuticals, nucleosides, steroids) When bonded to 300Å silica, it is an ideal phase for peptides, peptide mapping, and small hydrophilic proteins.	WAX	PEI, DEAE, Polyethylenimine, Diethylaminoethyl, Weak Base $-\text{Si}-\text{CH}_2\text{CH}_2\text{N}(\text{CH}_2\text{CH}_3)_2$ Ion exchange material. Weak anion exchangers (acidic) are most useful for analyzing acidic proteins and peptides.
C18	ODS, RP-18, LC18, Octadecyl $-\text{Si}-\text{C}_{18}\text{H}_{37}$ Classic Reversed Phase material. Most retentive for non-polar solutes. Excellent for ion-pairing chromatography. Wide applicability (e.g. nucleosides, nucleotides, steroids, pharmaceuticals, vitamins, fatty acids, environmental compounds). When bonded to 300Å silica, this phase is perfect for separating small hydrophilic peptides	WCX	CM, Carboxymethyl, Weak Acid $-\text{Si}-\text{CH}_2\text{COOH}$ Ion exchange material. Weak cation exchangers (basic) are most useful for analyzing basic proteins and peptides.
Phenyl, C ₆ H ₅	$-\text{Si}-\text{CH}_2\text{CH}_2\text{CH}_2$  Reversed Phase material. Unique selectivity. Useful for analyzing aromatic compounds. When bonded to 300Å silica, this phase is useful for HIC.		

SOLVENT MISCIBILITY TABLE

POLARITY CHART

SOLVENT	POLARITY INDEX	REFRACTIVE INDEX	UP (mm) D15	DOWN POINT (°C)	VELOCITY (g/min)	SOLUBILITY IN WATER (%)
Acetic Acid	6.2	1.372	230	118	1.26	100
Acetone	5.1	1.359	330	56	0.37	100
Acetonitrile	5.8	1.344	190	82	0.32	100
Benzene	2.7	1.501	280	80	0.65	0.18
Bulk Acetate	4.0	1.394	254	125	0.73	0.43
n-Butanol	3.9	1.399	215	118	2.98	7.61
Carbon tetrachloride	1.6	1.465	253	77	0.97	0.02
Chloroform	4.1	1.446	245	61	0.57	0.815
Cyclohexane	0.2	1.426	200	81	1.00	0.01
1,2-Dichloroethane	3.5	1.444	225	84	0.79	0.87
Dichloromethane	3.1	1.424	225	41	0.44	1.6
Dimethylformamide	6.3	1.431	268	155	0.92	100
Dimethyl sulfoxide	7.2	1.478	268	189	2.00	100
Dioxane	4.8	1.422	215	101	1.54	100
Ethyl Acetate	4.4	1.372	260	77	0.45	8.7
Ethanol	5.2	1.360	210	73	1.20	100
di-Ethyl Ether	2.8	1.353	220	35	0.32	6.89
Heptane	0.0	1.387	200	98	0.39	0.0003
Hexane	0.0	1.375	200	69	0.33	0.001
Methanol	5.1	1.329	205	65	0.60	100
Methoxy Ethyl Ether	2.5	1.369	210	55	0.27	4.8
Methyl Ethyl Ketone	4.7	1.379	328	80	0.45	24
Pentane	0.0	1.358	200	36	0.23	0.004
n-Propanol	4.0	1.384	210	97	2.27	100
iso-Propanol	3.9	1.377	210	82	2.30	100
iso-Propyl Ether	2.2	1.366	220	68	0.37	100
Tetrahydrofuran	4.0	1.407	215	65	0.55	100
Toluene	2.4	1.496	285	111	0.59	0.051
Trichloroethylene	1.0	1.477	273	97	0.57	0.11
Water	9.0	1.333	200	190	1.80	100
Xylene	2.5	1.500	290	139	0.61	0.018

Xylene	
Water	
Trichloroethylene	
Toluene	
Tetrahydrofuran	
Di-iso-Propyl Ether	
iso-Propanol	
n-Propanol	
pentane	
Methyl Ethyl Ketone	
Methyl-t-Butyl Ether	
Methanol	
Ethyl Acetate	
Hexane	
Dimethyl Sulfoxide	
Dimethylformamide	
Dichloromethane	
2-Dichloroethane	
Cyclohexane	
Chloroform	
Carbon tetrachloride	
Butanol	
Ethyl Acetate	
Heptane	
Octonitrile	
Octane	
Acetic Acid	

Immiscible
 Miscible

SYMBOLS TABLE
 1. Ethylene Oxide
 2. Methylene Chloride
 3. Methyl Sulfoxide
 4. tert-Butyl Methyl Ether
 5. 2-Butanone
 6. 2-Picoline

HPLC THEORY AND CALCULATIONS


HPLC CALCULATIONS

COLUMN EFFICIENCY: In general, N = Number of Theoretical Plates, a is a constant depending on method used, t_r = retention time of peak, and W = the peak width at a given peak height.

$$N = a \left(\frac{t_r}{W} \right)^2$$

Method	a
Peak Width $1/2$ Peak Height	5.54
Peak Width at 4.4% Peak Height (5σ method)	25
Tangent	16

The peak width at $1/2$ height is the most commonly used method for calculating HPLC column efficiency.

PEAK ASYMMETRY: $A_2 = B/A$ at 10% peak height.

CAPACITY FACTORS (RELATIVE RETENTION): The Capacity Factor, k' , of a sample component is a measure of the degree to which that component is retained by the column relative to an unretained component.

$$k' = (t_r - t_0) / t_0$$

Where t_r is the elution time of retained component, and t_0 is the elution time of the unretained sample.

SELECTIVITY: (α): $\alpha = k'_1 / k'_2$

RESOLUTION: R_s , defined as the amount of separation between two adjacent peaks, is given by

$$R_s = (1/4) (\alpha - 1) (N)^{1/2} (k'/1 + k')$$

where k' is the average value for the two peaks.

ADJUSTING FLOW RATE FOR DIFFERENT COLUMN IDS: When scaling up from analytical to preparative modes or when scaling down from analytical to microbore LC, it is often desirable to keep retention times constant. The flow rate can be adjusted so that the columns operate at the same linear velocity. When switching from a column with a radius (0.5 • I.D.) of r_1 to another with a radius of r_2 , the flow rate must be altered by a factor of X , where

$$X = (r_2 / r_1)^2$$

0.008

For example, when scaling up from a 250 x 4.6mm column to 250 x 10mm I.D. column, the flow rate must be increased by factor of 4.73 in the 10mm column to generate the same linear velocity as that of the 4.6mm I.D. column, as derived below

$$X = (5.0 / 2.3)^2 = 4.73$$

The general formula which will convert flow rate from any given column dimension to any other is as follows:


$$F_2 = F_1 \cdot (L_2 / L_1) \cdot (r_2 / r_1)^2$$

Where: L = length of the column, in mm
 r = radius of the column, in mm
 F = flow rate, in ml/min
 1 designates the first, or reference, column
 2 designates the second column

EFFECT OF DIFFERENT CONDITIONS ON SAMPLE RETENTION

Change in Separation	t_0	Effect on Retention Time	Band Spacing
Flow rate, F	$1/F$	$1/F$	None
Column volume, V_m	V_m	V_m	None
Increase in vol. % strong solvent	None	Decrease	Small change
New strong solvent	None	Changes	Changes
pH value	None	Changes	Changes
Column packing (e.g., cyano vs C18)	Little	Changes	Changes
Increase temperature	None	Decrease	Small change
New mobile phase additives	None	Changes	Changes

APPLICATIONS OF LIQUID CHROMATOGRAPHY

(From: D.L. Saunders, in *Chromatography*, 3rd ed, E. Heltmann, Ed., p. 81, Van Nostrand Reinhold, New York, 1975. With permission.)

OLUMN SELECTION GUIDE

les of Thumb